

21120.10120 SARS-CoV-2 Next Generation Sequencing Procedure

Copy of version 3.0 (approved and current)

Last Approval or Periodic Review Completed 02-Mar-2021

Controlled Copy ID 308478

Location client request

Next Periodic Review Needed On or Before 02-Mar-2022

Organization Eurofins - Viracor

Effective Date 09-Mar-2021

Comments for version 3.0

Removed (b) (4) from SOP because only an (b) (4) is plated for this assay.

Approval and Periodic Review Signatures

Type	Description	Date	Version	Performed By	Notes
Approval	(b) (6)	02-Mar-2021 11:06	3.0	(b) (6)	
				(b) (6)	
Approval	(b) (6) Approval	01-Mar-2021 21:03	3.0	(b) (6)	
				(b) (6)	
Approval	(b) (6) Approval	01-Mar-2021 17:02	3.0	(b) (6)	
				(b) (6)	
Approval	(b) (6)	18-Feb-2021 13:56	2.0	(b) (6)	
				(b) (6)	
Approval	(b) (6)	25-Jan-2021 13:16	2.0	(b) (6)	
				(b) (6)	

Approval (b) (6) 24-Jan-2021 20:21 2.0 (b) (6)

(b) (6)

Approval (b) (6) Approval 22-Jan-2021 18:47 2.0 (b) (6)

(b) (6)

Approval (b) (6) 12-Jan-2021 9:17 1.0 (b) (6)

(b) (6)

Approval (b) (6) Approval 09-Jan-2021 20:41 1.0 (b) (6)

(b) (6)

Approval (b) (6) Approval 05-Jan-2021 9:31 1.0 (b) (6)

(b) (6)

Version History

Version	Status	Type	Date Added	Date Effective	Date Retired
3.0	Approved and Current	Major revision	01-Mar-2021	09-Mar-2021	Indefinite
2.0	Retired	Major revision	22-Jan-2021	25-Jan-2021	09-Mar-2021
1.0	Retired	Initial version	04-Jan-2021	12-Jan-2021	25-Jan-2021

BioPharma Procedure for the SARS-CoV-2 Next-Generation Sequencing Assay

Table of Contents

INTENDED USE.....	2
TEST INFORMATION.....	2
ABBREVIATIONS.....	2
METHOD PRINCIPLE.....	2
SPECIMEN REQUIREMENTS.....	3
REAGENTS AND MATERIALS.....	4
CALIBRATORS AND STANDARDS.....	5
QUALITY CONTROL.....	6
EQUIPMENT AND SUPPLIES.....	8
PROCEDURE.....	8
CALCULATIONS.....	10
RESULT REPORTING AND REPEAT CRITERIA.....	10
EXPECTED VALUES AND REFERENCE RANGES.....	13
CLINICAL SIGNIFICANCE.....	13
PROCEDURE NOTES.....	13
LIMITATIONS OF THE METHOD.....	14
PERFORMANCE SPECIFICATIONS.....	14
SAFETY.....	14
RELATED DOCUMENTS.....	14
REFERENCES.....	15

INTENDED USE

The intended use of the assay is to provide SARS-CoV-2 nucleotide sequence determination for the S Gene by Next-Generation Sequencing.

TEST INFORMATION

Assay name: SARS-CoV-2 NGS

Date test initially placed into service: 25-Jan-2021

METHOD PRINCIPLE

Next Generation Sequencing

This assay uses conventional PCR chemistry for amplification of a SARS-CoV-2 target sequence from nucleic acids that have been extracted from nasopharyngeal swab specimens. The amplification step is followed by purification of the PCR products which is followed by visualization and concentration of the PCR products using the Agilent 2200/4200 TapeStation. Amplified samples are diluted, and libraries are prepared and indexed using the Illumina Nextera XT Library Prep Kit. Following a bead-based purification step, prepared libraries are quality checked on the Agilent 2200/4300 TapeStation. Libraries are then diluted, pooled, denatured, and loaded into a sequencing reagent cartridge for processing on the Illumina MiSeq instrument.

SPECIMEN REQUIREMENTS

Patient Preparation:

- No preparation is required.
- Specimen Collection and Transport

See SOP 21120.435 *Specimen Collection and Transport* for procedures for collecting the proper specimens and how to prepare the specimens for transportation to the laboratory.

Specimen Type and Handling

- NP Swab specimens received in saline from clients/sites/central lab are the acceptable specimen type for this assay.
- UTM/VTM specimens are proven stable for

(b) (4)

Extracted RNA should be stored frozen immediately following extraction if not preparing for RT-PCR within

(b) (4)

(b) (4)

REAGENTS AND MATERIALS

Description	Source	Part/Cat No	Storage/Expiration
Extraction			
EasyMAG Lysis Buffer	BioMerieux	280134	Ambient
EasyMAG Extraction Buffer 1		280130	
EasyMAG Extraction Buffer 2		280131	
EasyMAG Extraction Buffer 3	BioMerieux	280132	Unopened: Refrigerated Opened: Ambient
Magnetic Silica	BioMerieux	280133	Refrigerated
NGS Amplification			
Conventional RT-PCR primer mixes	(b) (4)		
(b) (4)			
D5000 ScreenTape	Agilent	5067-5588	Refrigerated
D5000 Sample Buffer	Agilent	5067-5589	Refrigerated
D5000 Ladder	Agilent	5067-5589	Refrigerated
Library Preparation			
Illumina Nextera XT Library Prep Kit - 24 samples (Box 1)	Illumina	FC-131-1024	-20°C
Illumina Nextera XT Library Prep Kit - 24 samples (Box 2)	Illumina	FC-131-1024	Refrigerated
(b) (4)			
Library Sequencing			
(b) (4)			
MiSeq v2 Reagent Kit – (b) (4)	Illumina	MS-(b) (4)	-20°C
(b) (4)			

QUALITY CONTROL

The quality control program for this test is established in accordance with SOP 21120.517 *Analytical Quality Control: Quality Control Procedures*.

Quality control samples are prepared and control ranges are established and maintained in accordance with SOP 21120.517 *Analytical Quality Control: Quality Control Procedures*.

Extraction Controls

Description	Source	Part/Cat No.	Storage / Expiration*
SARS-CoV-2 Whole Virus Positive Control	(b) (4)		

Description	Source	Part/Cat No.	Storage / Expiration*
-------------	--------	--------------	-----------------------

(b) (4)

*: See SOP 21120.380 *Expiration Dating of Laboratory Materials* for proper expiration dating assignment.

Control Procedure

SARS-CoV-2 Positive and Negative controls are included (b) (4)

Whole Virus Positive Control: This control has an approximate concentration of (b) (4)

Negative control:

(b) (4)

Quality Control Acceptance Criteria/Repeat Criteria

Thoroughly document all QC failures or repeat testing on the SARS-CoV-2 NGS Test Record 21120.10121.

Proficiency Testing

Refer to SOP 21120.384 *Proficiency Testing Program* for information on proficiency testing.

EQUIPMENT AND SUPPLIES

Description
Biological Safety Cabinet, Class II
Fixed-angle benchtop centrifuge (with rotor for 2mL tubes)
2mL non-skirted (b) (4) tubes
BioMerieux easyMAG
(b) (4) Thermal Cycler
(b) (4) PCR tubes or (b) (4) Optical 96-well Reaction Plate (standard plate)
(b) (4) caps or Optical Adhesive Covers
(b) (4) swing-bucket centrifuge
Agilent 2200 or 4200 TapeStation and associated components and software
(b) (4) Optical Tube Strips (8x)
(b) (4) Optical Tube Strip Caps (8x)
(b) (4) vortexer
PCR strip tube mini-centrifuge
(b) (4) 96-well MIDI plate magnetic stand
(b) (4) High-speed micro-plate shaker
96-well conical-bottom (b) (4) plate
(b) (4) Sequencer
Refrigerator capable of sustaining 2-8°C
Freezer capable of sustaining -15°C to -35°C
Freezer capable of sustaining -64°C to -90°C

PROCEDURE

Analytical

Nucleic Acid Extractions are performed using the one of the following SOPs:

- A. Refer to 21120.705 *NucliSense easyMAG Total Nucleic Acid Extration* for all instructions for nucleic acid isolations.

- 1.
- 2.
- 3.
- 4.
- 5.
- 6.
- 7.
- 8.
- 9.

(b) (4)

Conventional RT-PCR Procedure

Nucleic Acid Amplification is performed according to SOP 21120.461 *Real-Time PCR and RT-PCR Using* (b) (4) *Instruments*.

(b) (4)

(b) (4)

(b) (4)

(b) (4)

(b) (4)

(b) (4)

DNA Concentration Determination and Quality Check

Note: Purified product DNA should run as a single band or with very minor second band within an agarose gel or other liquid matrix. These matrices will reveal the presence of secondary products and contaminating DNAs and RNAs, but not proteins.

Remove D5000 ScreenTape(s), D5000 Sample Buffer and D5000 Ladder from storage at 2-8°C and allow reagents to equilibrate to room temperature for (b) (4).

(b) (4)

Refer to 21120.5983 *Agilent TapeStation Operation Maintenance and Calibration* for instructions on operation, maintenance, and calibration of the Agilent TapeStation instrument.

(b) (4)

10

Nextera XT Library Preparation (Tagmentation and Indexing)

9

(b) (4)

(b) (4)

Library Quality Check and Quantification

(b) (4)

Library Pooling

11

Library Denaturing and Dilution

MiSeq Prep and Loading

FDA-CBER-2022-1614-1790503

(b) (4)

Run Monitoring

Post-Run Wash

(b) (4)

Analysis

15

(b) (4)

(b) (4)

(b) (4)

(b) (4)

21120.369 Analytical Quality Control: Quality Control Procedures

21120.384 Proficiency Testing Program

21120.705 NucliSens easyMAG Total Nucleic Acid Extraction

21120.5983 Agilent Tapestation Operation Maintenance and Calibration

21120.443 *Thermal Cycler Preventative Maintenance Operation, Maintenance and Calibration of the Thermal Cycler* (b) (4)

21120.7375 *Illumina MiSeq Operation, Maintenance, and Calibration of the Illumina MiSeq / (b) (4) Sequencer*
Analytical Quality Control: Quality Control Procedures

21120.595 Specimen Processing Guide

21120.572 Verification and Qualification of Critical Laboratory Materials

21120.764 Preparation and QC of Oligonucleotides

21120.10121 SARS-CoV-2 NGS Test Record

21120.10181 SARS-CoV-2 Spike NGS Reagent Log

21120.6292 SARS-CoV-2 Sequencing Assay easyMAG Run Map

None